

THE UNIVERSITY INSTRUMENTATION CENTER IS PLEASED TO ANNOUNCE TWO NEW SERVICES



PIPETTE SERVICING and PRECISION WEIGHING

Service and repair of pipettes is performed on campus, under controlled conditions, using precision balances and OEM parts and following DIN 12650-6 and ISO 8655-6 standards.

Both services use the Mettler-Toledo MX5 balance with a measuring capability of 10^{-6} grams, repeatability of 0.8 μg , linearity of 1.0 μg , and an accuracy of plus or minus 2 micrograms. The balance is fully programmable with many time saving features.

Pipette Service includes:

- ✓ leakage testing,
- ✓ cleaning,
- ✓ maintenance,
- ✓ minor repairs,
- ✓ seal and O-ring replacement as needed using OEM parts,
- ✓ replacement of missing, broken or worn parts, and
 - ✓ extensive service reports (see back).
- ✓ Same day service is available for fast turn-around.

Pipette servicing is only \$10 per pipette, including preventive maintenance parts. Prior to service, the client must completely decontaminate all units of any radioactive or bio-hazardous materials.

Precision weighing services are available to students and the research community.

Please contact: Rob Cinq-Mars
University Instrumentation Center, 223 Parsons Hall
University of New Hampshire
Durham, New Hampshire 03824

Tel: 603-862-0174, Fax: 603-862-4104, email: Rob.Cinq-Mars@unh.edu, URL: www.unh.edu/uic